

A. Contents



1. Gastro⁺ Monitor
2. Operating Manual
3. D-Piece
4. Mouthpiece
5. (3) AA Batteries
6. (2) Stylus
7. Infection Control Guide
8. GastroCHART Tracking & Database Software CD & Manual
9. USB Cable
10. Cleaning Cloth
11. Carrying Case
12. Zippered Monitor Case
13. Carrying Strap
14. Hydrogen Breath Test Guide by Ledochowski

B. Monitor Overview



1. Display
2. On/Off Button
3. D-Piece Opening
4. D-Piece
5. Exhaust Port
6. Lug on D-Piece
7. Sounder
8. USB Connector
9. Mouthpiece
10. Battery Compartment
11. **Vent Hole** – DO NOT COVER WITH LABELS

**** NEW FEATURE v1.7****

A. Embedded software serial number





1. Touch gears  icon
2. Touch question  icon

Please Remember

- **Monitor must be calibrated before its first use.****
1. **Never** use or store alcohol based products including hand sanitizer around your monitor or tester.
 2. Change your D-piece every **30** days.
 3. Calibrate at least **monthly** to correct for sensor drift using **100 ppm** H₂ (Hydrogen) calibration gas.
 4. **Never** block the air/vent hole on the Gastro+ with labels etc.

C. Operating Instructions



1. Press and hold the ON/OFF button until the display becomes active. Release the button.
2. Insert the D-Piece into the monitor.
3. Place a mouthpiece into the D-Piece.
4. Start breath test | inhale and touch the  icon. (Patient should hold their breath for 15 seconds)
5. After the beep, exhale slowly into the mouthpiece, aiming to empty the lungs completely.
6. Monitor's display will show a rising ppm value.
7. During analyzation of the sample, an  hour glass will appear.
8. The highest ppm level will be displayed on the screen when the hour glass disappears and the  icon appears.
9. The reading will stay on the screen until you touch another icon.
10. To repeat the breath test, click the  icon and follow the step 5 above.
11. To switch the monitor off, press and hold the ON/OFF button for three seconds. The monitor will power off after 45 minutes of inactivity.

D. Indications



The Gastro⁺ can be used to investigate a variety of gastro-intestinal disorders including:

- Carbohydrate breakdown deficiency
- Carbohydrate malabsorption
- Lactose & Fructose intolerance
- Bacterial overgrowth
- Determination of time passage through the gut.

E. Infection Control



NEVER USE cleansing products containing *alcohol* or other cleaning agents or solvents around your monitor. coVita | Bedfont offers alcohol-free cleansing wipes which are specifically engineered for use on coVita | Bedfont monitors and accessories. Effective against 99.999% of Germs, Spores, Bacteria, and Fungi. Tested against H5N1 Avian Flu and H1N1 Swine Flu. Additionally, Bedfont offers a hand sanitizer which was designed to safely and gently cleanse your hands. Because it's engineered to be alcohol-free, it won't interfere with your monitor's sensor or accessories. Kills 99.99% of germs, bacteria, & fungi without drying your skin.